

Urine Drug Screen Cut Off Values

Amphetamine	1000 ng/ml
Barbiturate	200 ng/ml
Benzodiazepine	200 ng/ml
Cannabinoid	20 ng/ml
Cocaine	300 ng/ml
Fentanyl	5 ng/ml
Methadone	300 ng/ml
Opiates	300 ng/ml
Oxycodone	100 ng/ml
Propoxyphene	300 ng/ml

If the result is less than the cut-off value, the screen will be reported as negative.

If the result is greater than or equal to the cut off value, the screen will be reported as presumptive positive.

Specimens with presumptive positive results can be sent to our reference lab for confirmation by request only.